

A Complete Solution for Medical, Surgical & Healthcare Needs.



LOCKING PLATES	1
LOCKING SCREW	16
BONE PLATE	18
BONE SCREW	30
INTERLOCKING NAIL SYS	37
INTERLOCKING BOLT SYS	42
PROSTHESIS SYSTEM	44
SPINE SYSTEM	46





# SMALL COMPRESSION LOCKING PLATE 3.5MM

Catalog no.	011.155.03 to 011.155.16
Screw ø	3.5mm
Material	Stainless Steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



#### LCP SMALL T PLATE (OBLIQUE) 3.5MM - LEFT AND RIGHT

Catalog no.	011.151.02 to 011.151.10
Screw ø	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



# LCP SMALL T PLATE(RIGHT ANGLE) 3.5MM

Catalog no.	011.150.02 to 011.150.10
Screw ø	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



### LCP PROXIMAL HUMERUS PLATE 3.5MM (TYPE 1)

Catalog no.	011.156.04 to 011.156.12
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,6,8,10,12
<b>Proximal Screw</b>	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



#### LCP PROXIMAL HUMERUS PLATE 3.5MM (TYPE 2)

Catalog no.	011.1910.03 to 011.1910.10
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10
<b>Proximal Screw</b>	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



#### LCP PROXIMAL HUMERUS PLATE 3.5MM (TYPE 3)

Catalog no.	011.1911.03 to 011.1911.08
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L, Titanium
Holes	3,4,5,6,7,8
<b>Proximal Screw</b>	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)







#### LCP CLOVERLEAF PLATE 3.5MM

Catalog no.	011.154.03 to 011.154.03
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
<b>Proximal Screw</b>	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



## LCP RECONSTRUCTION PLATE 3.5MM ROUND HOLE

Catalog no.	011.153.03 to 011.153.18
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



## LCP RECONSTRUCTION PLATE 5.0MM (ROUND HOLE)

Catalog no.	011.1954.03 to 011.1954.18
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (Diff. 2.0mm)



### LCP RECONSTRUCTION PLATE 3.5MM (DYNAMIC HOLE)

011.1912.03 to 011.1912.16
3.5mm
Stainless steel 316L,316LVM, Titanium
3,4,5,6,7,8,9,10,11,12,13,14,15,16
Locking Screw 3.5mm
10mm to 40mm (Diff. 2.0mm)



## LCP RECONSTRUCTION PLATE 5.0MM (DYNAMIC HOLE)

Catalog no.	011.1913.03 to 011.1913.16
Screw ø	5.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (Diff. 2.0mm)



#### DISTAL TIBIA LOCKING PLATE LEFT & RIGHT

Catalog no.	011.157.04 to 011.157.12
Screw ø	5.0mm & 4.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4, 6, 8, 10, 12,
<b>Proximal Screw</b>	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking screw 5.0mm
Screw Length	
(Locked Screw 5	<b>5.0mm)</b> 10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5	<b>5.0mm)</b> 55mm to 90mm (Diff. 5.0mm)
· ·	







### MEDIAL DISTAL TIBIA LOCKING PLATE 3.5MM LEFT AND RIGHT

Catalog no.	011.1909.05 to 011.1909.13
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	5, 7, 9, 11, 13
<b>Proximal Screw</b>	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



#### DISTAL TIBIA ANATOMICAL LOCKING PLATE - LEFT & RIGHT

Catalog no.	011.1001.04 to 011.1001.14
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4, 6, 8, 10, 12, 14
<b>Proximal Screw</b>	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



# NARROW LOCKING COMPRESSION PLATE 5.0MM

Catalog no.	011.166.04 to 011.166.18
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Shaft Part	Locking Screw 5.0mm
Screw Length	
(Locked Screw §	<b>5.0mm)</b> 10mm to 50mm (Diff. 2.0mm)
(Locked Screw &	<b>5.0mm)</b> 55mm to 90mm (Diff. 5.0mm)



### BROAD LOCKING COMPRESSION PLATE 4.8MM

Catalog no.	011.167.04 to 011.167.20
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20
Shaft Part	Locking Screw 5.0mm
Screw Length	
(Locked Screw 5	<b>.0mm)</b> 10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5	<b>.0mm)</b> 55mm to 90mm (Diff. 5.0mm)



#### L BUTTRESS LOCKING PLATE LEFT & RIGHT

Catalog no.011.148.04 to 011.148.12Screw ø5.0mmMaterialStainless steel 316L,316LVM, TitaniumHoles4,5,6,7,8,9,10,11,12Proximal ScrewEast Screw LengthDistal ScrewLocking Screw 5.0mmScrew Length10mm to 50mm (Diff. 2.0mm)(Locked Screw 5.0mm)55mm to 90mm (Diff. 5.0mm)



#### **T PLATE LOCKING 5.0MM**

Catalog no.	011.162.04 to 011.162.12
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12
Proximal Screw	
Screw Length	
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locked Screw 5	<b>.0mm)</b> 10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5	<b>.0mm)</b> 55mm to 90mm (Diff. 5.0mm)







### T BUTTRESS LOCKING PLATE 5.0MM

Catalog no.	011.1262.04 to 011.1262.12
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12
Proximal Screw	
Screw Length	
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locked Screw 5	.0mm) 10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5	.0mm) 55mm to 90mm (Diff. 5.0mm)



#### DISTAL HUMERUS LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1925.03 to 011.1925.15
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,15
Shaft Part	Locking Screw 3.5mm
Screw Length	10.0mm to 60.0mm (diff. 2.0mm)



#### DISTAL HUMERUS WITH LATERAL SUPPORT LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1907.03 to 011.1907.15
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,15
Shaft Part	Locking Screw 3.5mm
Screw Length	10.0mm to 60.0mm (diff. 2.0mm)



### EXTRA-ARTICULAR DISTAL HUMERUS LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1941.04 to 011.1941.14
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,6,8,10,12,14
Shaft Part	Locking Screw 3.5mm
Screw Length	10.0mm to 60.0mm (diff. 2.0mm)



#### LATERAL TIBIAL LOCKING PLATE 5.0MM/6.5MM – LEFT & RIGHT

Catalog no.	011.161.03 to 011.161.12
Screw ø	5.0mm & 6.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12
<b>Proximal Screw</b>	Locking Screw 6.5mm
Screw Length	30mm to 95mm (Diff. 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locking Screw	5.0mm) 10mm to 50mm (Diff. 2.0mm)
(Locking Screw	<b>5.0mm)</b> 55mm to 90mm (Diff. 5.0mm)



#### LATERAL TIBIAL LOCKING PLATE 5.0MM – LEFT & RIGHT

Catalog no.	011.159.03 to 011.159.13
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
<b>Proximal Screw</b>	Cancellous Locking Screw 5.0mm
Screw Length	30mm to 80mm (Diff. 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locking Screw 5.0mm) 10mm to 50mm (Diff. 2.0mm)	
(Locking Screw	5.0mm) 55mm to 90mm (Diff. 5.0mm)







### DISTAL FEMUR LOCKING PLATE 5.0MM/6.5MM (LEFT AND RIGHT)

Catalog no.	011.163.03 to 011.163.13
Screw ø	5.0mm & 6.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
Proximal Screw	Cancellous Locking Screw 6.5mm
Screw Length	30mm to 95mm (Diff. 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locked Screw 5	.0mm) 10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5	.0mm) 55mm to 90mm (Diff. 5.0mm)



#### DISTAL FEMUR LOCKING PLATE LEFT & RIGHT

Catalog no.	011.164.03 to 011.164.13
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
<b>Proximal Screw</b>	Cancellous Locking Screw 5.0mm
Screw Length	30mm to 80mm (Diff. 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locking Screw 5.0mm) 10mm to 50mm (Diff. 2.0mm)	
(Locking Screw	5.0mm) 55mm to 90mm (Diff. 5.0mm)



## CONDYLAR BUTTERS LOCKING PLATE LEFT & RIGHT

Catalog no.	011.176.03 to 011.176.13
Screw ø	5.0mm & 6.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
<b>Proximal Screw</b>	Cancellous Locking Screw 6.5mm
Screw Length	30mm to 95mm (Diff. 5.0mm)
<b>Distal Screw</b>	Locking Screw 5.0mm
Screw Length	
(Locking Screw 5.0mm) 10mm to 50mm (Diff. 2.0mm)	
(Locking Screw	5.0mm) 55mm to 90mm (Diff. 5.0mm)



### DHS PLATE LOCKING LONG BARRAL 130<sup>o</sup> /135<sup>o</sup>

Catalog no.	011.165.02 to 011.165.12
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10,11,12
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)



#### DHS PLATE LOCKING SHORT BARRAL 130<sup>o</sup> /135<sup>o</sup>

Catalog no.	011.152.02 to 011.152.12
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10,11,12
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)



#### D.C.S PLATE LOCKING 90°

Catalog no.	011.19	937.03 to 011.1937.16
Screw ø	5.0mn	n & 6.5mm
Material	Stainle	ess steel 316L,316LVM, Titanium
Holes	3,4,5,6	6,7,8,9,10,11,12,13,14,15,16
Locking Screw 5	5.0mm	10mm to 50mm (Diff. 2.0mm)
		55mm to 90mm (Diff. 5.0mm)
Cancellous Locking		
Screw 6.5mm		30mm to 95mm (Diff. 5.0mm)







### DISTAL FEMUR LOCKING PLATE 5.0MM/6.5MM (LEFT AND RIGHT)

011.1928.04 to 011.1928.08
3.5mm
Stainless steel 316L,316LVM, Titanium
04, 06, 08
12mm to 40mm (Diff. 2.0mm)



### **PROXIMAL FEMUR LOCKING PLATE**

Catalog no.	011.1938.05 to 011.1938.13
Screw ø	5.0mm & 6.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	05, 07, 09, 11, 13
Locking Screw	
5.0mm	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
Cancellous	
Screw 6.5mm	30mm to 100mm (Diff. 5.0mm)



### PROXIMAL TIBIA LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1986.04 to 011.1986.14
Screw ø	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	04, 06, 08, 10, 12, 14
Locking Screw	
3.5mm	12mm to 40mm (Diff. 2.0mm)




#### LCP TOMIFIX MID HIGH T PLATE

Catalog no.	011.1252.03 & 011.1252.04
Screw ø	5.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	03, 04
Locking Screw	
5.0mm	18mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)



#### CALCANEUS LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1940.70mm & 011.1940.75mm
Screw ø	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	70mm, 75mm
Locking Screw	
3.5mm	12mm to 40mm (Diff. 2.0mm)



#### DISTAL TIBIA WITH TAB LOCKING PLATE - LEFT & RIGHT

Catalog no.	011.1967.05 to 011.1967.13
Screw ø	4.0 mm & 3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	05, 07, 09, 11, 13
Locking Screw	
3.5mm	10mm to 40mm (Diff. 2.0mm)
Locking Screw	
4.0mm	20mm to 70mm (Diff. 5.0mm)







### 2.4MM DISTAL RADIUS LOCKING PLATE STRAIGHT

Catalog no.	011.193.05 to 011.193.13
Screw ø	2.4mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	5, 7, 9, 11, 13
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



#### 2.4MM L DISTAL RADIUS LOCKING PLATE (HEAD 2HOLE) – LEFT & RIGHT

Catalog no.	011.191.3 & 011.191. 4
Screw ø	2.4mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	2+3, 2+4
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



#### 2.4MM L DISTAL RADIUS LOCKING PLATE (HEAD 3HOLE) – LEFT & RIGHT

Catalog no.	011.190. 3 & 011.190. 4
Screw ø	2.4mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3+3, 3+4
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



### 2.4MM L DISTAL RADIUS OBLIQUE LOCKING PLATE – LEFT & RIGHT

Catalog no.	011.194. 3 & 011.194. 4
Screw ø	2.4mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3+3, 3+4
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



#### 2.4MM T DISTAL RADIUS LOCKING PLATE (HEAD 3HOLE)

Catalog no.	011.192. 3 & to 011.192. 4
Screw ø	2.4mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3+3, 3+4
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



#### 2.4MM DISTAL RADIUS LOCKING PLATE (HEAD 5HOLE) – LEFT & RIGHT

Catalog no.	011.197. 3 & to 011.197. 5
Screw ø	2.4mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	5+3, 5+5
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)







### 2.4 MM BUTTRESS LOCKING PLATE (HEAD 5HOLE) – LEFT & RIGHT

Catalog no.	011.198
Screw ø	2.4mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	5+3
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



#### TWO-COLUMN DISTAL RADIUS PLATE LEFT & RIGHT

Catalog no.	011.1946.02 to 011.1946.04
Screw ø	2.4mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	02, 03, 04
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



## LOCKING SCREW

#### LCP SCREW DIA 2.4 MM

Catalog no	011.286.06 to 011.286.30
Thread Diameter	2.4mm
Use drills Diameter	2.0mm
Drive	hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	2.7mm
Length	06mm to 30mm (Diff. 2.0mm)

#### LCP SCREW DIA 2.7 MM

Catalog no	011.204.06 to 011.204.30
Thread Diameter	2.7mm
Use drills Diameter	2.0mm
Drive	hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	2.7mm
Length	06mm to 30mm (Diff. 2.0mm)



#### LCP SCREW DIA 3.5 MM

Catalog no	011.268.10 to 011.268.50
Thread Diameter	
Use drills Diameter	2.7mm
Drive	Hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	3.5mm
Length	10mm to 50mm (Diff. 2.0mm)



#### LCP SCREW DIA 4.0 MM

Catalog no	011.270.25 to 011.270.70
Thread Diameter	4.0mm
Use drills Diameter	3.2mm
Drive	Hex
Material	Stainless steel316L, 316LVM, Titanium
Size	4.0mm
Length	25mm to 70mm (Diff. 2.0mm)



## LOCKING SCREW





#### LCP SCREW DIA 5.0 MM

Catalog no	011.269.16 to 011.269.95
Thread Diameter	5.0mm
Use drills Diameter	4.0mm
Drive	Hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	5.0mm
Length	16mm to 50mm (Diff. 2.0mm)
	55mm to 95mm (Diff. 5.0mm)

### LCP SCREW CANCELLOUS DIA 5.0 MM

Catalog no	011.2930.30 to 011.2930.80
Thread Diameter	5.0mm
Use drills Diameter	4.0mm
Drive	Hex
Material	Stainless steel 316L,316LVM, Titanium
Size	5.0mm
Length	30mm to 50mm (Diff. 2.0mm)
	55mm to 80mm (Diff. 5.0mm)

#### LCP SCREW DIA 6.5 MM

Catalog no	011.271.35 to 011.271.95
Thread Diameter	6.5mm
Use drills Diameter	4.5mm
Drive	hex
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	35mm to 95mm (Diff. 5.0mm)



Desta and the



# RECONSTRUCTION PLATE FOR 2.7MM SCREW ROUND HOLE

Catalog no.	011.300.04 to 011.300.14
Screw ø	2.7mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 2.7mm
Screw Length	06mm to 30mm (Diff. 2.0mm)



#### DYNAMIC COMPRESSION PLATE FOR 2.7 MM SCREW

Catalog no.	011.320.04 to 011.320.14
Screw ø	2.7mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 2.7mm
Screw Length	06mm to 30mm (Diff. 2.0mm)



#### ONE THIRD TUBULAR PLATE FOR 3.5 MM SCREW

Catalog no.	011.321.04 to 011.321.14
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)







# RECONSTRUCTION PLATE FOR 3.5 MM SCREW ROUND HOLE

Catalog no.	011.322.04 to 011.322.18
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,13,14,15,16,17,18
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



## RECONSTRUCTION PLATE FOR 3.5 MM SCREW DC HOLE

Catalog no.	011.3221.04 to 011.3221.14
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,14
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



#### RECONSTRUCTION CURVED PLATE FOR 3.5MM SCREW

Catalog no.	011.301.04 to 011.301.16
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,08,09,10,11,12,13,14,15,16
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



# SMALL DCP PLATE FOR 3.5 MM SCREW

Catalog no.	011.325.04 to 011.325.14
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



## LC DCP PLATE SMALL FOR 3.5 MM SCREW

Catalog no.	011.333.04 to 011.333.14
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



### SMALL T PLATE FOR 3.5 MM SCREW

Catalog no.	011.326.03 to 011.326.10
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	03,04,05,06,07,08,09,10
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)







# SMALL T OBLIQUE PLATE FOR 3.5 MM SCREW

Catalog no.	011.327.03 to 011.327.10
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	03,04,05,06,07,08,09,10
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



100000

# ACRO PLATE FOR 3.5 MM SCREW LEFT & RIGHT

Catalog no.	011.3931.04 to 011.3931.06
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)

## RECONSTRUCTION PLATE FOR 4.5 MM SCREW ROUND HOLE

Catalog no.	011.323.04 to 011.323.18
Screw ø	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



# RECONSTRUCTION PLATE FOR 4.5 MM SCREW DC HOLE

Catalog no.	011.324.04 to 011.324.14
Screw ø	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14,
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



## SEMI TUBULAR PLATE FOR 4.5 MM SCREW

011.3225.04 to 011.3225.14
4.5mm
Stainless steel 316L,316LVM, Titanium
4,5,6,7,8,9,10,11,12,14,
Cortical Screw 4.5mm
16mm to 70mm (Diff. 2.0mm)



## NARROW DCP PLATE FOR 4.5 MM SCREW

Catalog no.	011.337.04 to 011.337.14
Screw ø	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)







# LC DCP PLATE NARROW FOR 4.5 MM SCREW

Catalog no.	011.334.04 to 011.334.14
Screw ø	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



## BROAD DCP PLATE FOR 4.5 MM SCREW

Catalog no.	011.338.04 to 011.338.18
Screw ø	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



## LC DCP PLATE BROAD FOR 4.5 MM SCREW

Catalog no.	011.358.06 to 011.358.14
Screw ø	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	6,7,8,9,10,11,12,14
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



### **CLOVERLEAF PLATE**

Catalog no.	011.328.03 to 011.328.12
Screw ø	
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12
Screw	
Screw Length	



#### **T PLATE FOR 4.5 MM SCREW**

Catalog no.	011.329.03 to 011.329.12
Screw ø	6.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12
<b>Proximal Screw</b>	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



## T BUTTRESS PLATE FOR 4.5 MM SCREW

Catalog no.	011.330.03 to 011.330.12
Screw ø	6.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12
<b>Proximal Screw</b>	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)







### L BUTTRESS PLATE FOR 4.5 MM SCREW LEFT & RIGHT

Catalog no.	011.332.03 to 011.332.12
Screw ø	6.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12
<b>Proximal Screw</b>	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)
<b>Distal Screw</b>	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



#### LATERAL TIBIAL PLATE (HOCKEY) FOR 4.5 MM SCREW LEFT & RIGHT

Catalog no.	011.335.04 to 011.335.14
Screw ø	6.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14
<b>Proximal Screw</b>	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



#### CONDYLAR BUTTRESS PLATE FOR 4.5 MM SCREW LEFT & RIGHT

Catalog no.	011.336.04 to 011.336.14
Screw ø	6.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14
<b>Proximal Screw</b>	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)
<b>Distal Screw</b>	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)





### DHS PLATE (DCP HOLES) SHORT BARREL

Catalog no.	011.339.04 to 011.339.14
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14
Degree	120 <sup>0</sup> , 125 <sup>0</sup> , 130 <sup>0</sup> , 135 <sup>0</sup> , 140 <sup>0</sup> , 145 <sup>0</sup>
<b>Proximal Screw</b>	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



#### DHS PLATE (ROUND HOLES) SHORT BARREL

Catalog no.	011.339.04R to 011.339.18R
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Degree	120 <sup>0</sup> , 125 <sup>0</sup> , 130 <sup>0</sup> , 135 <sup>0</sup> , 140 <sup>0</sup> , 145 <sup>0</sup>
<b>Proximal Screw</b>	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



#### DHS PLATE (DCP HOLES) LONG BARREL

Catalog no.	011.340.04 to 011.340.14
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14
Degree	120 <sup>0</sup> , 125 <sup>0</sup> , 130 <sup>0</sup> , 135 <sup>0</sup> , 140 <sup>0</sup> , 145 <sup>0</sup>
<b>Proximal Screw</b>	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)







### DHS PLATE (ROUND HOLES) LONG BARREL

Catalog no.	011.340.02R to 011.340.18R
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Degree	120 <sup>0</sup> , 125 <sup>0</sup> , 130 <sup>0</sup> , 135 <sup>0</sup> , 140 <sup>0</sup> , 145 <sup>0</sup>
<b>Proximal Screw</b>	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
<b>Distal Screw</b>	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



#### DHS PLATE (ROUND HOLES) PEDIATRICS

Catalog no.	011.331.04 to 011.331.10
Screw ø	
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10
Degree	125 <sup>0</sup> , 130 <sup>0</sup> , 135 <sup>0</sup>
<b>Proximal Screw</b>	Lag Screw
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



#### 95<sup>o</sup> DCS PLATE (DCP HOLES) TOTAL HOLE

Catalog no.	011.341.06 to 011.341.20
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	6,7,8,9,10,11,12,13,14,15,16,17,18,19,20
<b>Proximal Screw</b>	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
<b>Distal Screw</b>	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)
<b>Distal Screw</b>	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)



### DCS PLATE (ROUND HOLES) 95<sup>o</sup>

Catalog no.	011.342.06 to 011.342.20
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	6,7,8,9,10,11,12,13,14,15,16,17,18,19,20
<b>Proximal Screw</b>	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
<b>Distal Screw</b>	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)
<b>Distal Screw</b>	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)



## Y RECONSTRUCTION PLATE FOR 3.5mm SCREW

Catalog no.	011.302.01 to 011.302.10
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	1,2,3,4,5,6,7,8,9,10
Distal Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 4.0mm
Screw Length	10mm to 30mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 4.0mm
Screw Length	35mm to 70mm (Diff. 5.0mm)




### **CALCANEAL PLATE 3.5MM**

Catalog no.	011.310
Screw ø	4.0mm & 3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	Small, Medium, Large
Distal Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 4.0mm
Screw Length	10mm to 30mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 4.0mm







#### SCREW

Catalog no.	011.3944.04 to 10
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM ,Titanium
Holes	3,4,5,6,7,8,9,10,11,12
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)
Screw	Cortical Screw 4.5mm



011.311.04 to 011.311.14
Stainless steel 316L,316LVM, Titanium
4,5,6,7,8,9,10,11,12,13,14
Cortical Screw 4.5mm
16mm to 70mm (Diff. 2.0mm)
Cancellous Screw 6.5mm
30mm to 110mm (Diff. 5.0mm)

### TROCENTRIC PLATE FOR DHS

Catalog no.	011.3948
Screw ø	
Material	Stainless steel 316L,316LVM, Titanium
Holes	4 Hole
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)
<b>Proximal Screw</b>	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)

# 000

#### **EIGHT PLATE**

Catalog no.	011.366
Screw ø	
Material	Stainless steel 316L,316LVM, Titanium
Holes	4 Hole
Distal Screw	Cortical Screw 4.5mm





## CORTICAL SCREW (SELF TAPPING) 1.5MM

Catalog no.	011.49.06 to 011.49.30
Material	Stainless steel 316L,316LVM, Titanium
Size	1.5mm
Length	06mm to 30mm (Diff. 2.0mm)



### **CORTICAL SCREW**

#### 2.0MM

Catalog no.	011.49.06 to 011.49.30
Material	Stainless steel 316L,316LVM, Titanium
Size	1.5mm
Length	06mm to 30mm (Diff. 2.0mm)



### **CORTICAL SCREW**

#### 2.4MM

Catalog no.	011.485.06 to 011.485.30
Material	Stainless steel 316L,316LVM, Titanium
Size	2.4mm
Length	06mm to 30mm (Diff. 2.0mm)



### **CORTICAL SCREW 2.7MM**

Catalog no.	011.401.06 to 011.401.30
Material	Stainless steel 316L,316LVM, Titanium
Size	2.7mm
Length	06mm to 30mm (Diff. 2.0mm)







### **CORTICAL SCREW (14 TPI) 3.5MM**

Catalog no.	011.40
Material	Stainle
Size	3.5mm
Length	10mm

011.402.10 to 011.402.30 Stainless steel 316L,316LVM, Titanium 3.5mm 10mm to 60mm (Diff. 2.0mm)



### CORTICAL SCREW (14 TPI) SELF TAPPING 3.5MM

Catalog no.	011.01.402.10 to 011.01.402.60
Material	Stainless steel 316L,316LVM, Titanium
Size	3.5mm
Length	10mm to 60mm (Diff. 2.0mm)



# CORTICAL SCREW 3.5MM (20 TPI)

Catalog no.	011.403.10 to 011.403.10.60
Material	Stainless steel 316L,316LVM Titanium
Size	3.5mm
Length	10mm to 60mm (Diff. 2.0mm)



### CORTICAL SCREW (20 TPI) SELF TAPPING 3.5MM

Catalog no.	011.01.403.10 to 011.01.403.60
Material	Stainless steel 316L,316LVM, Titanium
Size	3.5mm
Length	10mm to 60mm (Diff. 2.0mm)



### CANCELLOUS SCREW (SHORT THREAD) 4.0MM

Catalog no.	011.404.10 to 011.404.70
Material	Stainless steel 316L,316LVM, Titanium
Size	4.0mm
Length	10mm to 30mm(Diff. 2.0mm)
	35mm to 70mm(Diff. 5.0mm)



### CANCELLOUS SCREW (FULL THREAD) 4.0MM

Catalog no.	011.405.10 to 011.405.70
Material	Stainless steel 316L,316LVM, Titanium
Size	4.0mm
Length	10mm to 30mm(Diff. 2.0mm)
	35mm to 70mm(Diff. 5.0mm)



### CANNULATED SCREW (SHORT THREAD) 4.0MM

Catalog no.	011.406.10 to 011.406.70
Material	Stainless steel 316L,316LVM, Titanium
Size	4.0mm
Length	10mm to 30mm(Diff. 2.0mm)
	35mm to 70mm(Diff. 5.0mm)



### CANNULATED SCREW (FULL THREAD) 4.0MM

Catalog no.	011.407.10 to 011.407.70
Material	Stainless steel 316L,316LVM, Titanium
Size	4.0mm
Length	10mm to 30mm(Diff. 2.0mm)
	35mm to 70mm(Diff. 5.0mm)







### **CORTICAL SCREW 4.5MM DIA**

Catalog no.	011.4
Vaterial	Stain
Size	4.5m
_ength	16mn

011.409.16 to 011.409.70 Stainless steel 316L,316LVM, Titanium 4.5mm 16mm to 70mm(Diff. 2.0mm)

## 

### CORTICAL SCREW (SELF TAPPING) 4.5MM

Catalog no.	011.01.409.12 to 011.01.409.70
Material	Stainless steel 316L,316LVM, Titanium
Size	4.5mm
Length	16mm to 70mm(Diff. 2.0mm)

### MALLEOLAR SCREW 4.5MM

Catalog no.	011.410.20 to 011.410.100
Material	Stainless steel 316L,316LVM, Titanium
Size	4.5mm
Length	20mm to 100mm(Diff. 5.0mm)



### **CANCELLOUS SCREW (16TH) 6.5MM**

Catalog no.	011.411.20 to 011.411.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)



### CANCELLOUS SCREW (32TH) 6.5MM

Catalog no.	011.412.20 to 011.412.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)



#### **CANCELLOUS SCREW (FULL TH)** 6.5MM

Catalog no.	011.413.25 to 011.413. 110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	25mm to 110mm(Diff. 5.0mm)

### **CANNULATED SCREW (16TH) 6.5MM**

Catalog no.	011.414.20 to 011.414.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)



### CANNULATED SCREW (32TH) 6.5MM

Catalog no.	011.415.25 to 011.415.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	25mm to 110mm(Diff. 5.0mm)



### **BONE SCREW**





#### **CANNULATED SCREW (FULL THREAD)** 6.5MM

Catalog no.	011.417.20
Material	Stainless ste
Size	6.5mm
Length	20mm to 10

to 011.417.110 eel 316L,316LVM, Titanium 0mm(Diff. 5.0mm)



#### **HARBAT SCREW**

Catalog no.	011.416.25 to 011.416.30
Material	Stainless steel 316L,316LVM, Titanium
Size	4.0mm
Length	25mm to 30mm(Diff. 2.0mm)



#### **DHS LAG SCREW**

Catalog no.	011.418.50 to 011.418.110
Material	Stainless steel 316L,316LVM, Titanium
Size	12.5mm
Length	50mm to 110mm(Diff. 5.0mm)



#### **DHS LAG SCREW (HEX)**

Catalog no.	011.01.418.50 to 011.01.418.110
Material	Stainless steel 316L,316LVM, Titanium
Size	12.5mm
Length	50mm to 110mm(Diff. 5.0mm)



#### DHS LAG SCREW PEDIATRICS

Catalog no.	011.463.40 to 011.463.105
Material	Stainless steel 316L,316LVM, Titanium
Size	
Length	40mm to 105mm(Diff. 5.0mm)



#### WASHER FOR CANCELLOUS SCREW

Catalog no.	011.4216
Material	Stainless steel 316L,316LVM, Titanium
Size	4.0mm, 7.0mm, 13mm
Length	Standard









Catalog no.011.545.20 to 011.545.48Locking Bolt øStainless steel 316L,316LVM ,TitaniumNail Diameter07mm, 08mm, 09mm, 10mm, 11mm, 12mmNail Length20cm to 46cm (Diff. 2.0cm)

#### **V NAIL FOR TIBIA**

Catalog no.	011.546.20 to 011.546.40
Locking Bolt ø	
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	07mm, 08mm, 09mm, 10mm, 11mm,
Nail Length	20cm to 40cm (Diff. 2.0cm)



#### SUPRA CONDYLAR NAIL

Catalog no.	011.548.15 to 011.548.30
Locking Bolt ø	
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	9mm,10mm,11mm,12mm
Nail Length	15cm to30cm(diff.5.0CM)
Proximal &	
Distal Bolt	
Bolt Length	



#### **RUSH NAIL**

Catalog no.	011.534.18 to 011.534.30
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	1.5mm, 2.0mm, 2.5mm, 3.0mm, 3.5mm,
	4.0mm, 4.5mm, 5.0mm, 6mm
Nail Length	18cm to 30cm (diff. 2.0cm)



#### SQUARE NAIL FOR ULNA

Catalog no.	011.547.2.0.18 to 011.547.2.0.30
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm,
	4.5mm, 5.0mm
Nail Length	18cm to 30cm (diff. 2.0cm)



#### SQUARE NAIL FOR RADIUS

Catalog no.	011.547.2.0.18 to 011.547.2.0.30
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm,
	4.5mm, 5.0mm
Nail Length	18cm to 30cm (diff. 2.0cm)



#### **HUMERUS INTERLOCKING NAIL**

Catalog no.	011.549.18 to 011.549.30
Locking Bolt ø	3.1mm
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	6.0mm, 7.0mm, 8.0mm
Nail Length	18cm to 30cm (diff. 2.0cm)
Proximal Bolt	Locking Bolt 3.1mm
Bolt Length	10mm to 40mm(diff. 2.0mm)
Distal Bolt	Locking Bolt 3.1mm
Bolt Length	10mm to 40mm(diff. 2.0mm)








#### **INTER LOCKING TIBIA NAIL (SOLID)**

Catalog no.	011.01.550.26 to 011.01.550.40
Locking Bolt ø	4.9mm
Material	Stainless steel 316L,316LVM ,Titanium
Nail Diameter	8.0mm, 9.0mm,10.0mm,11.0mm
Nail Length	26.0cm to 40.0cm(diff. 2.0cm)
Proximal &	
Distal Bolt	Locking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)
Proximal &	
Distal Bolt	Locking Bolt 4.9mm
Bolt Length	65mm to 80mm(diff. 5.0mm)



Catalog no.	011.550.26 to 011.550.40
Locking Bolt ø	4.9mm
Material	Stainless steel 316L,316LVM ,Titanium
Nail Diameter	8.0mm, 9.0mm,10.0mm,11.0mm
Nail Length	26.0cm to 40.0cm(diff. 2.0cm)
Proximal &	
Distal Bolt	Locking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)
Proximal &	
Distal Bolt	Locking Bolt 4.9mm
Bolt Length	65mm to 80mm(diff. 5.0mm)

#### **INTER LOCKING FEMUR NAIL**

Catalog no.	011.551.34. to 011.551.44
Locking Bolt ø	4.9mm
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	8.0mm, 9.0mm,10.0mm,11.0mm,12.0mm
Nail Length	34.0cm to 44.0cm (Diff. 2.0cm)
Proximal &	
Distal Bolt	Locking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)
Proximal &	
Distal Bolt	Locking Bolt 4.9mm
Bolt Length	65mm to 80mm(diff. 5.0mm)







#### **TFN NAIL SHORT**

Catalog no.	011.5935
Locking Bolt ø	8.0mm,6.4mm & 4.9mm
Material	Stainless steel 316L,316LVM Titanium
Nail Diameter	9.0mm, 10.0mm, 11.0mm, 12.0mm
Nail Length	18cm
Degree	130 <sup>o</sup> /135 <sup>o</sup>
Proximal Bolt	8.0mm
Bolt Length	40mm to 120mm(diff. 5.0mm)
Proximal Bolt	6.4mm
Bolt Length	35mm to 120mm(diff. 5.0mm)
Distal Bolt	Interlocking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)



#### **PFN NAIL SHORT**

Catalog no.	011.57
Locking Bolt ø	5.0mm, 3.7mm & 3.0mm
Material	Stainless steel 316L,316LVM Titanium
Nail Diameter	9.0mm, 10.0mm, 11.0mm, 12.0mm
Nail Length	25cm
Degree	130 <sup>o</sup> /135 <sup>o</sup>
Proximal Bolt	8.0mm
Bolt Length	40mm to 120mm(diff. 5.0mm)
Proximal Bolt	6.4mm
Bolt Length	35mm to 120mm(diff. 5.0mm)
Distal Bolt	Interlocking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)



#### PFN NAIL LONG 130<sup>o</sup> (LEFT & RIGHT)

011.57.09 to 011.57.12
3.0mm & 1.5mm
Stainless steel 316L,316LVM Titanium
9.0mm, 10.0mm, 11.0mm, 12.0mm
25cm
130º /135º
8.0mm
40mm to 120mm(diff. 5.0mm)
6.4mm
35mm to 120mm(diff. 5.0mm)
Interlocking Bolt 4.9mm
20mm to 60mm(diff. 2.0mm)
Left & Right







#### ELASTIC NAIL (FLEXIBLE NAIL)

Catalog no.	011.5938
Material	Stainless steel 316L,316LVM Titanium
Nail Diameter	2.0mm, 2.5mm, 3.0mm, 3.5mm,
	4.0mm, 4.5mm
Nail Length	44.0cm Standard

#### SIRUS NAIL (LEFT & RIGHT)



Catalog no.	011.559.34 to 011.559.44
Locking Bolt ø	6.5mm & 4.9mm
Material	Stainless steel 316L,316LVM Titanium
Nail Diameter	9.0mm,10.0mm,11.0mm,12.0mm
Nail Length	34cm to 44cm (diff. 2.0cm)
Proximal Bolt	Cannulated / Solid 6.5mm
Bolt Length	30mm to 110m (diff. 5.0mm)
Proximal &	
Distal Bolt	Locking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)
	65mm to 80mm(diff. 5.0mm)
Standard	Left & Right

#### **ENDER'S NAIL**

Catalog no.	011.553.18 to 011.553.46
Material	Stainless steel 316L,316LVM Titanium
Nail Diameter	3.0mm, 3.5mm, 4.0mm, 4.5mm,
Nail Length	18cm to 46cm (Diff. 1.0cm)



## **INTERLOCKING BOLT SYS**

#### **INTER LOCKING BOLT 3.1MM**

Catalog no.	011.633.18 to 011.633.50
Thread Diameter	3.1mm
Material	Stainless steel 316L,316LVM ,Titanium
Size	3.1mm
Length	18mm to 50mm (Diff. 2.0mm)



#### **INTER LOCKING BOLT 4.5MM**

Catalogue No.	011.01.652.18 to 011.01.652.90
Thread Diameter	4.5mm
Core Diameter	
Use drills Diameter	3.2mm
For passing hole	
Material	Stainless steel 316L,316LVM, Titanium
Size	4.5mm
Length	18mm to 90mm (Diff. 2.0mm)



#### **INTER LOCKING BOLT 3.9MM**

Catalog no.	011.653.18 to 011.653.80
Thread Diameter	3.9mm
Core Diameter	2.7mm
Use drills Diameter	3.2mm
For passing hole	4.0mm
Material	Stainless steel 316L,316LVM, Titanium
Size	3.9mm
Length	18mm to 50mm (Diff. 2.0mm)
	55mm to 80mm (Diff. 5.0mm)



#### **INTER LOCKING BOLT 4.9MM**

Catalog no.	011.652.18 to 011.652.90
Thread Diameter	4.9mm
Core Diameter	3.8mm
Use drills Diameter	4.0mm
For passing hole	5.1mm
Material	Stainless steel 316L,316LVM, Titanium
Size	4.9mm
Length	18mm to 60mm (Diff. 2.0mm)
	65mm to 90mm (Diff. 5.0mm)



### **INTERLOCKING BOLT SYS**





#### PFN BOLT CANNULATED 6.4MM

Catalog no.	011.65.30 to 011.65.120
Thread Diameter	
Core Diameter	
Use drills Diameter	
For passing hole	
Material	Stainless steel 316L,316LVM, Titanium
Size	6.4mm
Length	30mm to 120mm (Diff. 5.0mm)

#### **PFN BOLT CANCELLOUS 6.4MM**

Catalog no.	011.01.65.30 to 011.01.65.120
Thread Diameter	
Core Diameter	
Use drills Diameter	
For passing hole	
Material	Stainless steel 316L,316LVM, Titanium
Size	6.4mm
Length	30mm to 120mm (Diff. 5.0mm)

#### PFN BOLT CANNULATED 8.0MM

Catalog no.	011.66.45 to 120
Thread Diameter	
Core Diameter	
Use drills Diameter	
For passing hole	
Material	Stainless steel 316L,316LVM, Titanium
Size	8.0mm
Length	45mm to 120mm (Diff. 5.0mm)

#### PFN BOLT CANNULATED 8.0MM (HEX)

Catalog no.	011.01.66.45 to 011.01.66.120
Thread Diameter	
Core Diameter	
Use drills Diameter	
For passing hole	
Material	Stainless steel 316L,316LVM ,Titanium
Size	8.0mm
Length	45mm to 120mm (Diff. 5.0mm)









### **PROSTHESIS SYSTEM**

#### **THOMPSOM HIP PROSTHESIS**

Catlog no.	011.754.37 to 011.754.53
Material	Stainless steel 316L, 316LVM ,Titanium
Size	37, 39, 41, 43, 45, 47, 49, 51, 53
Length	Standard



#### **AUSTIN MOORE HIP PROSTHESIS**

Catlog no.	011.755.37 to 011.755.53
Material	Stainless steel 316L, 316LVM, Titanium
Size	37, 39, 41, 43, 45, 47, 49, 51, 53
Length	Standard



#### **BIPOLAR FIXED HEAD**

Catlog no.	011.78.37 to 011.78.53
Material	Stainless steel 316L, 316LVM, Titanium
Size	37, 39, 41, 43, 45, 47, 49, 51, 53
Length	Standard

#### **MODULAR BIPOLAR CUP**

Catlog no.	011.7115.37 to 011.7115.53
Material	Stainless steel 316L, 316LVM, Titanium
Size	37, 39, 41, 43, 45, 47, 49, 51, 53
Length	Small, Medium, Large

### **PROSTHESIS SYSTEM**



#### MODULAR BIPOLAR STEM

Catlog no.	-
Material	Stainless steel 316L, 316LVM, Titanium
Size	Small, Medium, Large
Length	Standard

#### SHOULDER PROSTHESIS

Catlog no.	011.751.35 to 011.751.51
Material	Stainless steel 316L, 316LVM, Titanium
Size	35, 37, 39, 41, 43, 45, 47, 49, 51
Length	Standard

#### AUSTIN MOORE HIP PROSTHESIS LONG STEM

Catlog no.	011.01.755.35 to 011.01.755.51
Material	Stainless steel 316L, 316LVM, Titanium
Size	35, 37, 39, 41, 43, 45, 47, 49, 51
Length	Standard







#### MONO AXIAL SCREW

011.826.20 to 011.826.50
Titanium
4.5mm, 5.5mm, 6.5mm, 7.5mm
20mm to 50mm(Diff. 5.0mm)



#### **POLY AXIAL SCREW**

Catalog no.	011.828.20 to 011.828.50
Material	Titanium
Dia	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)



#### MONO AXIAL SCREW (D/T)

Catalog no.	011.833.20 to 011.833.50
Material	Titanium
Dia	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)



#### POLY AXIAL SCREW (D/T)

Catalog no.	011.834.20 to 011.834.50
Material	Titanium
Dia	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)







#### MONO AXIAL REDUCTION SCREW

Catlog no.	011.827.20 to 011.827.50
Material	Titanium
Size	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)



#### POLY AXIAL REDUCTION SCREW

Catlog no.	011.829.20 to 011.829.50
Material	Titanium
Size	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)



#### MONO AXIAL REDUCTION (D/T)

Catlog no.	011.813.20 to 011.813.50
Material	Titanium
Size	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)



#### POLY AXIAL REDUCTION (D/T)

Catlog no.	011.895.20 to 011.895.50
Material	Titanium
Size	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)



#### 3.5 MM DIA - MONOAXIAL SCREW

011.830.12 to 011.830.40
Titanium
3.5mm, 4.2mm, 5.0mm, 6.0mm
12mm to 20mm(Diff. 2.0mm)
25mm to 40mm(Diff. 5.0mm)



#### 3.5MM DIA – POLYAXIAL SCREW

Catalog no.	011.831.12 to 011.831.40
Material	Titanium
Dia	3.5mm, 4.2mm
Length	12mm to 20mm(Diff. 2.0mm)
	25mm to 40mm(Diff. 5.0mm)



#### LAMINAR HOOK ANGLE BLADE

Catalog no.	011.835.01
Material	Titanium
Dia	Standard
Length	Standard



#### LAMINAR HOOK NARROW BLADE

Catalog no.	011.831.12 to 011.831.40
Material	Titanium
Dia	3.5mm, 4.2mm
Length	12mm to 20mm(Diff. 2.0mm)
	25mm to 40mm(Diff. 5.0mm)








#### **PEDICAL HOOK**

Catlog no.	011.835.03
Material	Titanium
Size	Standard
Length	Standard



#### LEMINAR HOOK WIDE BLADE

Catlog no.	011.835.04
Material	Titanium
Size	Standard
Length	Standard



#### STAPLE

Catlog no.	011.829.01
Material	Titanium
Size	Standard
Length	Standard

#### CROSS LINK CONNECTORE 5MM X 3MM ROD

Catlog no.	011.838.01
Material	Titanium
Size	Standard
Length	Standard



#### 3.5mm Dia - CONNECTING ROD PRE BEND

Catalog no. Material Dia Length 011.823.3.5 Titanium 3.5mm Standard

#### **CONNECTING ROD**

Catalog no.	011.832
Material	Titanium
Dia	3.0mm, 3.5mm, 5.0mm, 5.5mm
Length	50mm to 450mm



#### 3.5mm Dia - POSTERIOR OCCIPITAL PLATE WITH ROD

Catalog no.	011.854.3.5.100 to 011.854.3.5.300
Material	Titanium
Dia	3.5mm
Length	100mm, 200mm, 250mm, 300mm



#### **POSTERIOR OCCIPITAL SCREW**

Catalog no.	011.8133.06 to 011.8133.20
Material	Titanium
Dia	3.5mm
Length	06mm to 20mm (Diff. 2.0mm)







# ANTERIOR CERVICAL CAGE (SQUARE TYPE)

Catlog no.	011.844.4 to 8.14 to 11.844.4 to 8.16
Material	Titanium
Size	4.0mm, 5.0mm, 06mm, 7.0mm, 8.0mm
Length	12mmx14mm, 12mmx16mm



#### **EXPANDABLE CAGE**

Catlog no.	011.846.20 to 011.846.75
Material	Titanium
Size	10mm, 12mm, 14mm, 16mm, 18mm, 20mm
Length	20mm to 75mm (Diff. 5.0mm)



#### ECO CAGE

Catlog no.	011.845.15 to 011.845.100
Material	Titanium
Size	10mm, 12mm, 14mm, 16mm, 18mm, 20mm
Length	15mm to 100mm (Diff. 5.0mm)



#### **PLIF CAGE**

Catlog no.	011.847.20 to 011.847.30
Material	Titanium
Size	7.0mm, 8.0mm, 9.0mm, 10.0mm, 11.0mm
Length	20mm, 25mm, 30mm



#### TLIF CAGE (BANANA TYPE)

Catalog no.	011.856.7 to 011.856.11
Material	Titanium
Dia	7.0mm, 8.0mm, 9.0mm, 10.0mm, 11.0mm
Length	40, 80, 110



#### PLIF CAGE (BULLET TYPE)

 Catalog no.
 011.857.20 to 011.857.30

 Material
 Titanium

 Dia
 8.0mm, 10.0mm, 12.0mm

 Length
 20mm, 25mm, 30mm



#### ANTERIOR CERVICAL PLATE (Z TYPE)

Catalog no.	011.842.25 to 011.842.70
Material	Titanium
Dia	Standard
Length	25.0mm,27.5mm,30.0mm,32.5mm,
	35.0mm,37.5mm,40.0mm,42.5mm,
	45.0mm,47.5mm,50.0mm,55.0mm,
	60.0mm,65.0mm,70.0mm



#### BONE SCREW FOR ACP 4mm

#### TITANIUM

Catalog no.	011.843.10 to 011.843.24
Material	Titanium
Dia	4.0mm
Length	10mm to 24mm (Diff. 2.0mm)



#### **OCCIPITAL LOOPS**

Catlog no.	011.855
Material	Titanium
Size	4.0mm
Length	



## 2.4MM LCP SET

ITEM DESCRIPTION	QTY
DISTAL RADIUS LOCKING PLATE 2.4, STRAIGHT, 5, 6 HOLES	4
L DISTAL RADIUS LOCKING PLATE 2.4	-
LEFT ANGLED, SHAFT 3, 4 HOLES, HEAD 2 HOLES	4
RIGHT ANGLED, SHAFT 3, 4 HOLES, HEAD 2 HOLES	4
LEFT ANGLED, SHAFT 3, 4 HOLES, HEAD 3 HOLES	4
RIGHT ANGLED, SHAFT 3, 4 HOLES, HEAD 3 HOLES	4
OBLIQUE, LEFT ANGLED, SHAFT 3, 4 HOLES, HEAD 3 HOLES	4
OBLIQUE, RIGHT ANGLED, SHAFT 3, 4 HOLES, HEAD 3 HOLES	4
T DISTAL RADIUS LOCKING PLATE 2.4, SHAFT 3, 4 HOLES, HEAD 3 HOLES	4
DISTAL RADIUS LOCKING PLATE 2.4, LEFT, SHAFT 3, 5 HOLES, HEAD 5 HOLES	4
DISTAL RADIUS LOCKING PLATE 2.4, RIGHT, SHAFT 3, 5 HOLES, HEAD 5 HOLES	4
BUTTRESS LOCKING PLATE 2.4, LEFT, SHAFT 3 HOLES, HEAD 5 HOLES	2
BUTTRESS LOCKING PLATE 2.4, RIGHT, SHAFT 3 HOLES, HEAD 5 HOLES	2
LCP VOLAR 2.4MM TWO-COLUMN DISTAL RADIUS PLATE LEFT 02, 03, 04HOLE	6
LCP VOLAR 2.4MM TWO-COLUMN DISTAL RADIUS PLATE RIGHT 02, 03, 04HOLE	6
2.4MM LCP SCREW 06MM TO 30MM (10 EACH)	260
2.4MM CORTICAL SCREW 06MM TO 30MM (10 EACH)	260
LCP DRILL SLEEVE 2.0MM	2
DRILL BIT 2.MMX150MM	1
2.4MM SELF HOLDING SCREW DRIVER	1
2.4MM SCREW DRIVER	1

# **3.5MM RECONSTRUCTION LCP SET**

ITEM DESCRIPTION	QTY
LCP RECONSTRUCTION PLATE FOR 3.5MM SCREW ROUND HOLE 04HOLE	4
LCP RECONSTRUCTION PLATE FOR 3.5MM SCREW ROUND HOLE 05HOLE	4
LCP RECONSTRUCTION PLATE FOR 3.5MM SCREW ROUND HOLE 06HOLE	4
LCP RECONSTRUCTION PLATE FOR 3.5MM SCREW ROUND HOLE 07HOLE	4
LCP RECONSTRUCTION PLATE FOR 3.5MM SCREW ROUND HOLE 08HOLE	4
LCP RECONSTRUCTION PLATE FOR 3.5MM SCREW ROUND HOLE 09HOLE	4
LCP RECONSTRUCTION PLATE FOR 3.5MM SCREW ROUND HOLE 10HOLE	4
LCP RECONSTRUCTION PLATE FOR 3.5MM SCREW ROUND HOLE 11HOLE	4
LCP RECONSTRUCTION PLATE FOR 3.5MM SCREW ROUND HOLE 12HOLE	4
LCP RECONSTRUCTION PLATE FOR 3.5MM SCREW ROUND HOLE 13HOLE	4
LCP RECONSTRUCTION PLATE FOR 3.5MM SCREW ROUND HOLE 14HOLE	4
3.5MM LOCKING SCREW 10MM TO 30MM (5 EACH)	110
3.5MM CORTICAL SCREW SELF TAPPIING (20 TPI) 10MM TO 30MM (5 EACH)	110
DRILL SLEEVE 2.5MM	2
DRILL BIT 2.5MMX200MM	1
TAP 3.5MM (20 TPI)	1
DEPTH GUAGE 3.5MM	1
SCREW DRIVER 3.5MM	1
BOX FOR RECONSTRUCTION LCP SET	1

# **5.0MM RECONSTRUCTION LCP SET**

ITEM DESCRIPTION	QTY
LCP RECONSTRUCTION PLATE FOR 5.0MM SCREW ROUND HOLE 04HOLE	4
LCP RECONSTRUCTION PLATE FOR 5.0MM SCREW ROUND HOLE 05HOLE	4
LCP RECONSTRUCTION PLATE FOR 5.0MM SCREW ROUND HOLE 06HOLE	4
LCP RECONSTRUCTION PLATE FOR 5.0MM SCREW ROUND HOLE 07HOLE	4
LCP RECONSTRUCTION PLATE FOR 5.0MM SCREW ROUND HOLE 08HOLE	4
LCP RECONSTRUCTION PLATE FOR 5.0MM SCREW ROUND HOLE 09HOLE	4
LCP RECONSTRUCTION PLATE FOR 5.0MM SCREW ROUND HOLE 10HOLE	4
LCP RECONSTRUCTION PLATE FOR 5.0MM SCREW ROUND HOLE 11HOLE	4
LCP RECONSTRUCTION PLATE FOR 5.0MM SCREW ROUND HOLE 12HOLE	4
LCP RECONSTRUCTION PLATE FOR 5.0MM SCREW ROUND HOLE 13HOLE	4
LCP RECONSTRUCTION PLATE FOR 5.0MM SCREW ROUND HOLE 14HOLE	4
5.0MM LOCKING SCREW 20MM TO 50MM (5 EACH)	160
4.5MM CORTICAL SCREW SELF TAPPIING 20MM TO 50MM (5 EACH)	160
DRILL SLEEVE 4.0MM	2
DRILL BIT 4.0MMX200MM	1
DRILL BIT 3.2MMX200MM FOR CORTICAL SCREW	1
TAP 4.5MM	1
DEPTH GUAGE 4.5MM	1
SCREW DRIVER 4.5MM	1
BOX FOR RECONSTRUCTION LCP SET	1

## **BROAD LCP SET**

ITEM DESCRIPTION	QTY
LCP BROAD COMPRESSION PLATE FOR 5.0MM SCREW 04HOLE	2
LCP BROAD COMPRESSION PLATE FOR 5.0MM SCREW 05HOLE	2
LCP BROAD COMPRESSION PLATE FOR 5.0MM SCREW 06HOLE	2
LCP BROAD COMPRESSION PLATE FOR 5.0MM SCREW 07HOLE	2
LCP BROAD COMPRESSION PLATE FOR 5.0MM SCREW 08HOLE	2
LCP BROAD COMPRESSION PLATE FOR 5.0MM SCREW 09HOLE	2
LCP BROAD COMPRESSION PLATE FOR 5.0MM SCREW 10HOLE	2
LCP BROAD COMPRESSION PLATE FOR 5.0MM SCREW 11HOLE	2
LCP BROAD COMPRESSION PLATE FOR 5.0MM SCREW 12HOLE	2
LCP BROAD COMPRESSION PLATE FOR 5.0MM SCREW 13HOLE	2
LCP BROAD COMPRESSION PLATE FOR 5.0MM SCREW 14HOLE	2
LCP BROAD COMPRESSION PLATE FOR 5.0MM SCREW 16HOLE	2
5.0MM LOCKING SCREW 20MM TO 50MM 55/60/65/70 (5 each)	200
4.5MM CRTICAL SCREW SELF TAPPIING 20MM TO 56MM 55/60/65/70 (5 each)	200
DRILL SLEEVE 4.0MM	2
DRILL BIT 4.0MMX200MM	1
DRILL BIT 3.2MMX200MM FOR CORTICAL SCREW	1
TAP 4.5MM	1
DEPTH GUAGE 4.5MM	1
SCREW DRIVER 4.5MM	1
BOX FOR BROAD LCP SET	1

### **CEMENTED BIPOLAR SET - MODULAR**

ITEM DESCRIPTION	QTY
STEM WITH ONE CENTRALISER (small, medium, large, XL, XXL)	5
CUP 37MM TO 53MM S/M/L - 1 EACH	27
CEMENT RESTRICTOR (STERILE)	1
BIPOLAR FIX 37MM TO 53MM - 1 EACH	9
Instrument Set	1

### **CLAVICLE LCP SET**

ITEM DESCRIPTION	QTY
CLAVICLE LOCKING PLATE - LEFT 6 8 10 HOLE	6
CLAVICLE LOCKING PLATE - RIGHT 6 8 10 HOLE	6
3.5MM LOCKING SCREW 10MM TO 40MM (5 EACH)	160
3.5MM CORTICAL SCREW SELF TAPPIING (20 TPI) 10MM TO 40MM (5 EACH)	160
DRILL SLEEVE 2.5MM	2
DRILL BIT 2.5MMX200MM	1
TAP 3.5MM (20 TPI)	1
DEPTH GUAGE 3.5MM	1
SCREW DRIVER 3.5MM	1
BOX FOR CLAVICLE LCP SET	1

## **CLOVERLEAF PLATE LCP SET**

ITEM DESCRIPTION	QTY
LCP CLOVERLEAF PLATE, LENGTH 03HOLE TO 13HOLE (1 each)	22
3.5MM LOCKING SCREW 10MM TO 40MM (5 each)	160
3.5MM CRTICAL SCREW SELF TAPPIING (20 TPI) 10MM TO 40MM (5 each)	160
DRILL SLEEVE 2.5MM	2
DRILL BIT 2.5MMX200MM	1
TAP 3.5MM (20 TPI)	1
DEPTH GUAGE 3.5MM	1
SCREW DRIVER 3.5MM	1
BOX FOR CLOVERLEAF LCP SET	1

## DHS DCS SET

ITEM DESCRIPTION	QTY
DHS PLATE LCP 130D & 135D LONG	
04HOLE, 05HOLE, 06HOLE, 07HOLE, 08HOLE, 09HOLE, 10HOLE, 12HOLE (2 EACH)	16
DCS PLATE LCP 95D	
06HOLE, 07HOLE, 08HOLE, 09HOLE, 10HOLE, 12HOLE (1 EACH)	6
DHS LAG SCREW 50MM TO 115MM (1 EACH)	14
5.0MM LOCKING SCREW 20MM TO 50MM (5 EACH)	80
4.5MM CORTICAL SCREW SELF TAPPIING 20MM TO 50MM (5 EACH)	80
GUIDE PIN 2.5 MM WITH THREAD TIP	2
DHS ANGLE GUIDE - VARIABLE	1
T-HANDLE WITH QUICK COUPLING	1
SCREW DRIVER 3.5MM TIP	1
DHS TRIPPLE REAMER	1
DCS TRIPLE REAMER	1
WRENCH FOR DHS SCREWS	1
DHS IMPACTOR	1
DHS TAP	1
CENTERING SLEEVE FOR DHS TAP	1
CENTERING SLEEVE FOR DHS WRENCH	1
COUPLING SCREW FOR INSERTING DHS SCREW	1
COMPRESSOR	1
DCS ANGLE GUIDE	1
DHS INSTRUMENT SET	1

# DISTAL FEMUR LCP SET

ITEM DESCRIPTION	QTY
LCP CLOVERLEAF PLATE, LENGTH 03HOLE TO 13HOLE (1 each)	22
3.5MM LOCKING SCREW 10MM TO 40MM (5 each)	160
3.5MM CRTICAL SCREW SELF TAPPIING (20 TPI) 10MM TO 40MM (5 each)	160
DRILL SLEEVE 2.5MM	2
DRILL BIT 2.5MMX200MM	1
TAP 3.5MM (20 TPI)	1
DEPTH GUAGE 3.5MM	1
SCREW DRIVER 3.5MM	1
BOX FOR CLOVERLEAF LCP SET	1

# DISTAL HUMERUS LCP SET

ITEM DESCRIPTION	QTY
LCP DISTAL HUMERUS PLATE FOR 3.5MM/2.7MM SCREW	12
LEFT 03H / 05H / 07H / 09H / 11H / 13H (1 EACH)	
LCP DISTAL HUMERUS PLATE FOR 3.5MM/2.7MM SCREW	12
RIGHT 03H / 05H / 07H / 09H / 11H / 13H (1 EACH)	
LCP DISTAL HUMERUS WITH LATERAL SUPPORT FOR 3.5MM/2.7MM SCREW	12
LEFT 03H / 05H / 07H / 09H / 11H / 13H (1 EACH)	
LCP DISTAL HUMERUS WITH LATERAL SUPPORT FOR 3.5MM/2.7MM SCREW	12
RIGHT 03H / 05H / 07H / 09H / 11H / 13H (1 EACH)	
LCP DISTAL HUMERUS MEDIAL FOR 3.5MM/2.7MM SCREW	12
LEFT 03H / 05H / 07H / 09H / 11H / 13H (1 EACH)	
LCP DISTAL HUMERUS MEDIAL FOR 3.5MM/2.7 SCREW	12
RIGHT 03H / 05H / 07H / 09H / 11H / 13H (1 EACH)	
2.7MM LOCKING SCREW 10MM TO 40MM (5 EACH)	160
3.5MM LOCKING SCREW 10MM TO 40MM (5 EACH)	160
3.5MM CORTICAL SCREW SELF TAPPIING 10MM TO 40MM (5 EACH)	160
DRILL SLEEVE 2.5MM	2
DRILL SLEEVE 3.2MM	2
DRILL BIT 4.0MMX200MM	1
TAP 3.5MM	1
DEPTH GUAGE 3.5MM	1
SCREW DRIVER 3.5MM	1
BOX FOR DISTAL HUMERUS LCP SET	1

### DISTAL TIBIA LCP SET

ITEM DESCRIPTION	QTY
LCP DISTAL TIBIA PLATE FOR 5.0MM/4.0MM SCREW LEFT 04H / 06H / 08H / 12H	8
LCP DISTAL TIBIA PLATE FOR 5.0MM/4.0MM SCREW RIGHT 04H / 06H / 08H / 12H	8
5.0MM LOCKING SCREW 20MM TO 50MM (5 EACH)	160
4.5MM CORTICAL SCREW SELF TAPPIING 20MM TO 50MM (5 EACH)	160
4.0MM LOCKING SCREW 20MM TO 70MM (4 EACH)	88
4.0MM CANCELLOUS SCREW FULL TH 30MM TO 80MM (4 EACH)	88
DRILL SLEEVE 4.0MM	2
DRILL SLEEVE 3.2MM	2
DRILL BIT 4.0MMX200MM	1
DRILL BIT 3.2MMX200MM	1
TAP 4.5MM	1
TAP 4.0MM C.C	1
DEPTH GUAGE 4.5MM	1
SCREW DRIVER 3.5MM	1
SCREW DRIVER 4.5MM	1
BOX FOR DISTAL TIBIA LCP SET	1

# INTERLOCKING HUMERUS SET

ITEM DESCRIPTION	QTY
INTERLOCKING HUMRUS NAIL 6MM, 7MM, 8MM X 18CM TO 30CM	63
3.1mm interlocking bolt 10mm to 30mm (2 Each)	66
INSTRUMENT BOX	1
BONE AWL	1
САР	1
CONICAL BOLT	2
DEPTH GUAGE	1
DRILL BIT 2.5MMX250MM	1
DRILL SLEEVE 2.5MM	1
EXTRACTOR ROAD	1
EXTRACTOR HANDLE	1
GUIDE WIRE BOL TIP	2
GUIDE WIRE PLAIN	2
HAMMER	1
PROTECTION SLEEVE	1
SCREW DRIVER 3.5MM	1
SPENNAR	1
HUMERUS ZIG	1
TISSUE PROTECTOR	1
TOMEY BAR	1
TORCAR	1

# L BUTTRESS LCP SET

ITEM DESCRIPTION	QTY
LCP L BUTTRESS PLATE FOR 5.0MM SCREW LEFT 03HOLE TO 12HOLE	20
LCP L BUTTRESS PLATE FOR 5.0MM SCREW RIGHT 03HOLE TO 12HOLE	20
5.0MM LOCKING SCREW 20MM TO 50MM (5 each)	160
4.5MM CORTICAL SCREW SELF TAPPIING 20MM TO 50MM (5 each)	160
DRILL SLEEVE 4.0MM	2
DRILL BIT 4.0MMX200MM	1
DRILL BIT 3.2MMX200MM FOR CORTICAL SCREW	1
TAP 4.5MM	1
DEPTH GUAGE 4.5MM	1
SCREW DRIVER 4.5MM	1
BOX FOR L BUTTRESS LCP SET	1

# LATERAL TIBIAL LCP SET

ITEM DESCRIPTION	QTY
LCP LATERAL TIBIAL PLATE (HOCKEY) (AO) LEFT 04HOLE TO 13HOLE (1each)	10
LCP LATERAL TIBIAL PLATE (HOCKEY) (AO) RIGHT 04HOLE TO 13HOLE (1each)	10
5.0MM LOCKING SCREW 20MM TO 50MM (5 each)	160
5.0MM LOCKING SCREW 55MM TO 90MM (4 each)	64
4.5MM CRTICAL SCREW SELF TAPPIING 20MM TO 58MM (5 each)	200
5.0MM LOCKING CANSLOUSH SCREW 30MM TO 80MM (4 each)	88
DRILL SLEEVE 4.0MM	2
DRILL BIT 4.0MMX200MM	1
TAP 4.5MM	1
DEPTH GUAGE 4.5MM	1
SCREW DRIVER 4.5MM	1
BOX FOR LATERAL TIBIAL LCP SET	

# MEDIAL DISTAL TIBIA LCP SET

ITEM DESCRIPTION	QTY
LCP MEDIL DISTAL TIBIA PLATE LEFT 05H / 07H / 09H / 11H / 13H (1 each)	10
LCP MEDIL DISTAL TIBIA PLATE RIGHT 05H / 07H / 09H / 11H / 13H (1 each)	10
3.5MM LOCKING SCREW 10MM TO 40MM (5 each)	160
3.5MM CRTICAL SCREW SELF TAPPIING 10MM TO 40MM (5 each)	160
4.0MM LOCKING SCREW 20MM TO 70MM (4 each)	88
4.0MM CANCELLOUS SCREW FULL TH 30MM TO 80MM (4 each)	88
DRILL SLEEVE 2.5MM	2
DRILL SLEEVE 3.2MM	2
DRILL BIT 2.5MMX200MM	1
DRILL BIT 3.2MMX200MM	1
TAP 3.5MM	1
TAP 4.0MM C.C	1
DEPTH GUAGE 3.5MM	1
SCREW DRIVER 3.5MM	1
BOX FOR MEDIL DISTAL TIBIA LCP SET	1

# NARROW LCP SET

ITEM DESCRIPTION	QTY
LCP NARROW COMPRESSION PLATE FOR 5.0MM SCREW 04HOLE	2
LCP NARROW COMPRESSION PLATE FOR 5.0MM SCREW 05HOLE	2
LCP NARROW COMPRESSION PLATE FOR 5.0MM SCREW 06HOLE	2
LCP NARROW COMPRESSION PLATE FOR 5.0MM SCREW 07HOLE	2
LCP NARROW COMPRESSION PLATE FOR 5.0MM SCREW 08HOLE	2
LCP NARROW COMPRESSION PLATE FOR 5.0MM SCREW 09HOLE	2
LCP NARROW COMPRESSION PLATE FOR 5.0MM SCREW 10HOLE	2
LCP NARROW COMPRESSION PLATE FOR 5.0MM SCREW 11HOLE	2
LCP NARROW COMPRESSION PLATE FOR 5.0MM SCREW 12HOLE	2
LCP NARROW COMPRESSION PLATE FOR 5.0MM SCREW 13HOLE	2
LCP NARROW COMPRESSION PLATE FOR 5.0MM SCREW 14HOLE	2
5.0MM LOCKING SCREW 20MM TO 50MM 55/60/65/70 (5 each)	200
4.5MM CORTICAL SCREW SELF TAPPIING 20MM TO 50MM (5 each)	200
DRILL SLEEVE 4.0MM	2
DRILL BIT 4.0MMX200MM	1
DRILL BIT 3.2MMX200MM FOR CORTICAL SCREW	
TAP 4.5MM	
DEPTH GUAGE 4.5MM	
SCREW DRIVER 4.5MM	
BOX FOR NARROW LCP SET	

## **PFN INSTRUMENT SET**

ITEM DESCRIPTION	QTY
PFN NAIL SHORT 130D & 135D,09MM,10MM,11MM,12MM X 25CM (1EACH)	24
PFN NAIL LONG 130D (LEFT & RIGHT) 09MM,10MM,11MM X 34CM TO 44CM	108
6.4MM PFN BOLT 30MM TO 120MM (2 EACH)	114
8.0MM PFN BOLT 40MM TO 120MM (2 EACH)	102
4.9MM INTERLOCKING BOLT 28MM TO 60MM (2 EACH)	102
BOX FOR PFN INSTRUMENT SET	1
PFN C.C.AWL	1
PFN SOLAPUR SLEEVE	1
PFN CAP FOR CONICAL BOLT	1
PFN CONICAL BOLT - 2	2
PFN DRILL BIT 4.0MMX250MM	1
PFN DRILL SLEEVE 4.0MM	1
PFN EXTRACTOR HANDLE	1
PFN EXTRACTOR ROD	1
PFN GUIDE WIRE 1.8MMX400MM	2
PFN GUIDE WIRE 2.0MMX400MM	2
PFN GUIDE WIRE SLEEVE 1.8MM	1
PFN GUIDE WIRE SLEEVE 2.0MM	1
PFN HAMMER	1
PFN ZIG 130D	1
PFN ZIG 135D	1

## **PFN INSTRUMENT SET**

ITEM DESCRIPTION	QTY
PFN PROTECTION SLEEVE (4.9MM)	1
PFN PROTECTION SLEEVE (6.4MM)	1
PFN PROTECTION SLEEVE (8.0MM)	1
PFN REMAR 6.4MM	1
PFN REMAR 8.0MM	1
PFN SCREW DRIVER 4.5MM	1
PFN SCREW DRIVER C.C 6.4MM	1
PFN SPENAR FOR CONICAL BOLT	1
PFN SPENAR FOR EXTRACTOR HANDLE	1
PFN T HANDLE QUICKCOUPLING FOR REMAR	1
PFN T WRENCH FOR 8.0MM SCREW	1
PFN TOMAY BAR	1
PFN TROCAR 6.4MM	1
PFN TROCAR 8.0MM	1

# PROXIMAL FEMUR LCP SET

ITEM DESCRIPTION	QTY
PROXIMAL FEMUR LOCKING PLATE 05,07,09,11,13 (1 each)	10
5.0MM LOCKING SCREW 20MM TO 50MM (5 each)	160
5.0MM LOCKING SCREW 55MM TO 90MM (4 each)	64
4.5MM CORTICAL SCREW SELF TAPPIING 20MM TO 58MM (5 each)	200
6.5MM CANCELLOUS SCREW FULL TH 30MM TO 95MM (2 each)	56
DRILL SLEEVE 4.0MM	2
DRILL BIT 4.0MMX200MM	1
DRILL BIT C.C 4.5MMX200MM	1
TAP 4.5MM	1
TAP 6.5MM C.C	1
DEPTH GUAGE 4.5MM	1
SCREW DRIVER 4.5MM	1
BOX FOR PROXIMAL FEMUR LCP SET	1

# PROXIMAL HUMERUS (NEW) LCP SET

ITEM DESCRIPTION	QTY
LCP PROXIMAL HUMERUS PLATE AO, LENGTH 03HOLE TO 10HOLE (1 each)	32
3.5MM LOCKING SCREW 10MM TO 40MM (5 each)	320
3.5MM CORTICAL SCREW SELF TAPPIING (20 TPI) 10MM TO 40MM (5 each)	320
4.0MM LOCKING SCREW 20MM TO 70MM (4 each)	176
4.0MM CANCELLOUS SCREW FULL TH 30MM TO 80MM (4 each)	176
DRILL SLEEVE 2.5MM	2
DRILL SLEEVE 3.2MM	2
DRILL BIT 2.5MMX200MM	1
DRILL BIT 3.2MMX200MM	1
TAP 3.5MM (20 TPI)	1
TAP 4.0MM C.C	1
DEPTH GUAGE 3.5MM	1
SCREW DRIVER 3.5MM	1
BOX FOR PROXIMAL HUMERUS (NEW) LCP SET	1

### PROXIMAL TIBIA LCP SET

ITEM DESCRIPTION	QTY
PROXIMAL TIBIA LOCKING PLATE - LEFT 04,06,08,10,12 HOLE (1 each)	10
PROXIMAL TIBIA LOCKING PLATE - RIGHT 04,06,08,10,12 HOLE (1 each)	10
3.5MM LOCKING SCREW 10MM TO 40MM (5 each)	160
3.5MM CORTICAL SCREW SELF TAPPIING 10MM TO 40MM (5 each)	160
4.0MM LOCKING SCREW 20MM TO 70MM (4 each)	88
4.0MM CANCELLOUS SCREW FULL TH 30MM TO 80MM (4 each)	88
DRILL SLEEVE 2.5MM	2
DRILL SLEEVE 3.2MM	2
DRILL BIT 2.5MMX200MM	1
DRILL BIT 3.2MMX200MM	1
TAP 3.5MM	1
TAP 4.0MM C.C	1
DEPTH GUAGE 3.5MM	1
SCREW DRIVER 3.5MM	1
BOX FOR PROXIMAL DISTAL LCP SET	1

# SMALL COMPRESSION SET

ITEM DESCRIPTION	QTY
LCP SMALL COMPRESSION PLATE FOR 3.5 MM SCREW, LENGTH 04HOLE	4
LCP SMALL COMPRESSION PLATE FOR 3.5 MM SCREW, LENGTH 05HOLE	4
LCP SMALL COMPRESSION PLATE FOR 3.5 MM SCREW, LENGTH 06HOLE	4
LCP SMALL COMPRESSION PLATE FOR 3.5 MM SCREW, LENGTH 07HOLE	4
LCP SMALL COMPRESSION PLATE FOR 3.5 MM SCREW, LENGTH 08HOLE	4
LCP SMALL COMPRESSION PLATE FOR 3.5 MM SCREW, LENGTH 09HOLE	4
LCP SMALL COMPRESSION PLATE FOR 3.5 MM SCREW, LENGTH 10HOLE	4
LCP SMALL COMPRESSION PLATE FOR 3.5 MM SCREW, LENGTH 11HOLE	4
LCP SMALL COMPRESSION PLATE FOR 3.5 MM SCREW, LENGTH 12HOLE	4
3.5MM LOCKING SCREW 10MM TO 40MM (5 each)	320
3.5MM CORTICAL SCREW SELF TAPPIING (20 TPI) 10MM TO 40MM (5 each)	320
DRILL SLEEVE 2.5MM	2
DRILL BIT 2.5MMX200MM	1
TAP 3.5MM (20 TPI)	1
DEPTH GUAGE 3.5MM	1
SCREW DRIVER 3.5MM	1
BOX FOR SMALL LCP SET	1

# SMALL T & T OBLIQUE LCP SET

ITEM DESCRIPTION	QTY
LCP SMALL T PLATE RIGHT ANGLE, LENGTH 03HOLE TO 10HOLE (1 each)	32
LCP SMALL T PLATE OBLIQUE LEFT, LENGTH 03HOLE TO 10HOLE (1 each)	32
LCP SMALL T PLATE OBLIQUE RIGHT, LENGTH 03HOLE TO 10HOLE (1 each)	32
3.5MM LOCKING SCREW 10MM TO 30MM (5 each)	220
3.5MM CORTICAL SCREW SELF TAPPIING (20 TPI) 10MM TO 30MM (5 each)	220
DRILL SLEEVE 2.5MM	2
DRILL BIT 2.5MMX200MM	1
TAP 3.5MM (20 TPI)	1
DEPTH GUAGE 3.5MM	1
SCREW DRIVER 3.5MM	1
BOX FOR SMALL T & T OBLIQUE LCP SET	1

### SPINE INSTRUMENT SET

ITEM DESCRIPTION	QTY
MONOAXIAL SCREW (DLX) 5.5MM X 30,35,40,45MM (4 each)	16
MONOAXIAL SCREW (DLX) 6.5MM X 30,35,40,45MM (4 each)	16
POLYAXIAL SCREW (DLX) 5.5MM X 30,35,40,45MM (4 each)	16
POLYAXIAL SCREW (DLX) 6.5MM X 30,35,40,45MM (4 each)	16
CONNECTIG ROD 5.5 MM	24
CONNECTING ROD 3.0MM	20
CONNECTOR	4
STAPLES	4
PROBE FLAT	1
PROBE FLAT CVD	1
TROCAR 3MM	1
TROCAR 4MM	1
BONE AWL STR	1
SOUND 2MM	1
SOUND 3MM	1
MONOAXIAL SCREW DRIVER	1
POLYAXIAL SCREW DRIVER	1
REDUCTION POLYAXIAL SCREW DRIVER	1
VICE GRIP	1
ROD HOLDING 5MM	1
ROD HOLDING 3MM	1
ROD PUSHER OPEN TYPE	1
ROD PUSHER CLOSED TYPE	1
COMPRESSOR CVD	1
DISTRACTOR STR	1

# T BUTTRESS LCP SET

ITEM DESCRIPTION	QTY
LCP T BUTTRESS PLATE FOR 5.0MM SCREW 03HOLE TO 12HOLE	20
5.0MM LOCKING SCREW 20MM TO 50MM (5 each)	160
4.5MM CORTICAL SCREW SELF TAPPIING 20MM TO 50MM (5 each)	160
DRILL SLEEVE 4.0MM	2
DRILL BIT 4.0MMX200MM	1
DRILL BIT 3.2MMX200MM FOR CORTICAL SCREW	1
TAP 4.5MM	1
DEPTH GUAGE 4.5MM	1
SCREW DRIVER 4.5MM	1
BOX FOR T BUTTRESS LCP SET	1

# T PLATE LCP SET

ITEM DESCRIPTION	QTY
LCP T PLATE FOR 5.0MM SCREW 03HOLE TO 12HOLE	20
5.0MM LOCKING SCREW 20MM TO 50MM (5 each)	160
4.5MM CORTICAL SCREW SELF TAPPIING 20MM TO 50MM (5 each)	160
DRILL SLEEVE 4.0MM	2
DRILL BIT 4.0MMX200MM	1
DRILL BIT 3.2MMX200MM FOR CORTICAL SCREW	1
TAP 4.5MM	1
DEPTH GUAGE 4.5MM	1
SCREW DRIVER 4.5MM	1
BOX FOR T PLATE LCP SET	1

# **TFN INSTRUMENT SET**

ITEM DESCRIPTION	QTY
TFN NAIL SHORT 130D & 135D, 09MM,10MM,11MM,12MM X 18CM (1each)	24
6.4MM PFN BOLT 30MM TO 120MM (2 each)	114
8.0MM PFN BOLT 40MM TO 120MM (2 each)	102
4.9MM INTERLOCKING BOLT 28MM TO 60MM (2 each)	102
BOX FOR PFN INSTRUMENT SET	1
PFN C.C.AWL	1
PFN CAP FOR CONICAL BOLT	1
PFN CONICAL BOLT - 2	2
PFN DRILL BIT 4.0MMX250MM	1
PFN DRILL SLEEVE 4.0MM	1
PFN EXTRACTOR HANDLE	1
PFN EXTRACTOR ROD	1
PFN GUIDE WIRE 1.8MMX400MM	2
PFN GUIDE WIRE 2.0MMX400MM	2
PFN GUIDE WIRE SLEEVE 1.8MM	1
PFN GUIDE WIRE SLEEVE 2.0MM	1
PFN HAMMER	1
PFN ZIG 130D	1
PFN ZIG 135D	1
PFN PROTECTION SLEEVE (4.9MM)	1
PFN PROTECTION SLEEVE (6.4MM)	1
PFN PROTECTION SLEEVE (8.0MM)	1

# **TFN INSTRUMENT SET**

ITEM DESCRIPTION	QTY
PFN REAMAR 6.4MM	1
PFN REAMAR 8.0MM	1
PFN SCREW DRIVER 4.5MM	1
PFN SCREW DRIVER C.C 6.4MM	1
PFN SPENAR FOR CONICAL BOLT	1
PFN SPENAR FOR EXTRACTOR HANDLE	1
PFN T HANDLE QUICKCOUPLING FOR REMAR	1
PFN T WRENCH FOR 8.0MM SCREW	1
PFN TOMAY BAR	1
PFN TROCAR 6.4MM	1
PFN TROCAR 8.0MM	1

### **TIBIA & FEMUR NAILING SET**

ITEM DESCRIPTION	QTY
INTER LOCKING (GUNDRILL) TIBIA NAIL 9MM,10MM,11MM X 28CM TO 38CM (1 each)	90
INTER LOCKING (GUNDRILL) FEMUR NAIL 9MM,10MM,11MM,12MM X 34CM TO 44CM (1 each)	120
4.9MM INTERLOCKING BOLT 22MM TO 60MM, 65MM, 70MM (2 each)	220
INSTRUMENT SET WITH COMPLETE INSTRUMENT (Distal Zig)	1
AWL	1
САР	1
CONICAL BOLT- 2	2
DEPTH GUAGE	1
DRILL BIT 4.0MMX200MM	1
DRILL BIT 3.2MMX250MM	1
DRILL SLEEVE 3.2MM	1
DRILL SLEEVE 4.0MM	
EXTRACTOR ROD	
EXTRACTOR HANDLE	
GUIDE WIRE BOL TIP	
GUIDE WIRE PLAIN	
HAMMER	
PROTECTION SLEEVE	
SCREW DRIVER 4.5MM	
SPENNAR	
TIBIA & FEMUR ZIG	
DISTAL FOR FEMUR	
DISTAL FOR TIBIA	
TISSUE PROTECTOR	
TOMEY BAR	
TORCAR	

#### **UNCEMENTED BIPOLAR SET**

ITEM DESCRIPTION	QTY
ATI- BIOACTIV VG STEM (12/14) TAPER	
SIZE-EX SMALL/SMALL+/MEDIUM/MEDIUM+/LARGE	6
BIPOLAR SHELL 2 PREFIX LINER WITH LOCK RING 37 TO 53	
SIZE-SHORT/STANDARD/LONG( (-5/0/+3)	9
FEMORAL HEAD Ø CC 28 MM (12/14) CONE	
SIZE-SHORT/STANDARD/LONG/EX LONG( (-5/0/+3+7)	4
FEMORAL HEAD Ø CC 28 MM (12/14) CONE	
SIZE-SHORT/STANDARD/LONG/EX LONG( (-5/0/+3+7)	2
INSTRUMENT SET	1

#### **UNCEMENTED THR SET**

ITEM DESCRIPTION	QTY
ATI- BIOACTIV VG STEM (12/14) TAPER H.A. N PORUS COATED	12
SIZE-EX SMALL/SMALL+/MEDIUM/MEDIUM+/LARGE	
BIOACTIV ACETABULAR TI UR SHELL	14
SIZE-42/44/46/48/50/52/54/56/58/60/62	
ID 22/28 MM URDR ACETABULAR CUP	12
SIZE-32/36/40/44/48 MM	
FEMORAL HEAD Ø CC 28 MM (12/14) CONE	8
SIZE-SHORT/STANDARD/LONG/EX LONG( (-5/0/+3+7)	
FEMORAL HEAD Ø CC 22 MM (12/14) CONE SIZE-STANDARD/LONG (TWO EACH)	4
UR SCREW TI SIZE-16/20/25/30/35/40/45/50/55/60 MM	24
INSTRUMENT SET	

**CONSOPHARMAPLUS** is an emerging manufacturer-exporter based in New Delhi, capital city of India. The company is founded by highly experienced professionals from Medical Devices and Healthcare industry.

The company has a vision to provide high quality Orthopedic solutions and other healthcare products at a very competitive cost so that the each and every class of society gets quality treatment and care. We ensure our customers and users of our products for the best value of their money.

Our main products include Orthopedic Implants, Orthopedic Instruments, Orthopedic Sets, Surgical Equipments and Rehabilitation Aids etc. The company has production facilities in Gujarat, India which is well known for its industrial expertise. We have completely automated and state-of-art manufactuing unit for orthopedic implants and instruments which assures quality control as per international standards.

Apart from Orthopedic Solutions, we offer our customers a complete solution for their Medical, Surgical and Healthcare needs. These product line include Hospital Furniture, Medical Disposables, Laboratory Supplies, Autoclaves and sterilizers, Suction Units, other medical devices and healthcare equipments.

We welcome dealers, doctors, healthcare professionals to experience all time best services and co-operation. Come, let us work together for serving the community called "Humanity" – a step towards making our world more healthy and beautiful !

CONSOPHARMA PLUS RZ-C 50, Nihal Vihar, Nangloi, New Delhi - 110041 India. Call : +91.935050 4074, 987300 4074 Email: consopharmaplus@gmail.com