



CONSOPHARMA PLUS®



*A Complete Solution for  
Medical, Surgical & Healthcare Needs.*





## Plain Hospital Bed - Deluxe

- Perforated CRCA sheet top
- S.S. bows of 32 mm dia.
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated



## Semi Fowler Bed with Metal Panel

- Two section perforated CRCA sheet top
- Back rest operated by screw with folding handles.
- Detachable S.S. bows of 32 mm dia. with metal panel.
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated



## Fowler Bed

- Four section perforated CRCA sheet top
- Back rest and Knee rest operated by screw with folding handles.
- Detachable S.S. bows of 32 mm dia. with metal panel, 4 flush I.V. rod locations
- SS swing type railing (optional)
- 125 m dia heavy duty wheels; 2 with brakes
- Pre-treated and epoxy powder coated

## ICU Bed - Electric

- Ventilated trapeze type bed surface
- All four movements operated by wired remote
- Backrest movement 80°
- Trendelenburg movement 12°
- Collapsible railings
- Buffers at four corners 150 mm dia,  
4 heavy duty wheels, 2 with stopper



## Hospital Bed with Baby Cradle

- Ventilated trapeze type bed surface
- All four movements operated by wired remote
- Backrest movement 80°
- Trendelenburg movement 12°
- Collapsible railings
- Buffers at four corners 150 mm dia,  
4 heavy duty wheels, 2 with stopper



## Pediatric Bed

- Complete with full length dropside safety railings
- Head and leg bows of equal height
- Pre-treated and epoxy powder coated





## Baby Cradle

- M.S. tubular pipe.
- Stand fitted with PVC stumps.
- Wire mesh hanging baby cradle.
- Pre-treated and Epoxy Powder coated



## Labour Table (Telescopic)

- Three section Stainless Steel Top.
- Back rest on ratchet; head low/up by screw mechanism.
- Middle Section with 'U' notch with bucket/ bowl.
- Sliding foot end. Railings on three sides.
- Available in full stainless steel / stainless steel top & epoxy powder coated base
- Lithotomy Bars Gutter type / Leather Strip type



## Emergency Trolley

- Height adjustable trolley with X-Ray permeable detachable top and full length X-Ray cassette holder
- Height adjustment, trendelenburg & reverse trendelenburg positions by screw and lever mechanism
- 150 mm dia heavy duty wheels, 2 with brakes
- Provision for fixing I.V. rod
- Swing away type S.S. railings on both side, rubber buffers on all corners
- Pre-treated and epoxy powder coated

## Stretcher on Trolley

- Detachable Stainless steel top fixed on CRCA tube frame mounted on 150 mm dia heavy duty wheels
- Provision for fixing I.V. rod
- Pre-treated and epoxy powder coated



## Examination Couch

- 4 Cabinets; 1 with lock
- Drawer with lock
- Sliding stepping tool & Writing pad
- 100 mm thick form mattress.
- Frame made of CRCA sheets and rectangular tube, CRCA sheet cabinet
- Pre treated & Epoxy Powder coated



## Kick Bucket and Bowl Stand

- Two layer Stainless Steel bowl stand
- 2 Stainless steel bowls of dia 345 mm each
- Stand available in stainless steel/ epoxy powder coated





## Foot Step - Double

- M.S. Tubular legs fitted with PVC stump.
- M.S. top fixed with rubber matting and aluminum beadings.
- Pre-treated and Epoxy Powder Coated.



## Revolving Stool

- Stainless steel top
- Completely made in SS
- Top base size 300 mm dia
- Height adjustable from 510 mm to 710 mm



## Emergency Trolley

- 6 Drawer Modular system upper two drawers with medicine containers
- 6 Detachable plastic bins
- Cages for oxygen and nitrous cylinders
- Top tray for monitor and pulse oxymeter.
- System mounted on 100 mm dia. wheels, with pair of front brakes.

## Bedside Table - Deluxe

- CRCA sheet main cabinet and a drawer.
- Mounted on 50 mm dia. twin roller wheels, front wheels with brake.
- Stainless steel tops with 3 sides raised and front turned down.
- Pre-treated and epoxy powder coated.



## Over Bed Table/Food Trolley

- Stainless steel top
- Height adjustable by lever gear mechanism
- Mounted on 50 mm dia wheels
- Pre-treated and epoxy powder coated



## Mayo's Trolley

- Trolley Mounted on 50 mm dia. wheels
- Available in full stainless steel / Powder coated frame with stainless steel tray
- Tray Height adjustable from 860 mm to 1300 mm







## Instrument Trolley

- Stainless steel bottom and top shelves with railing on top shelf
- Trolley Mounted on 100 mm dia. wheels
- Available in full stainless steel / Powder coated frame with stainless steel shelves and railings



## Dressing Trolley

- 2 stainless steel shelves with railing on top
- Supplied with Stainless steel bowl and bucket
- Available in full stainless steel/ Powder coated frames with stainless steel shelves and railings

## Drip Stand

- Available with two/ four hooks
- Made of full stainless steel/ epoxy powder coated stand & stainless steel saline rod.
- Five winged SS/PP base mounted on 5cm Castor.
- Height adjustable from Min.:1370mm to Max.:2435mm



## Bed Side Screen (3 Panels)

- Tubular frame work mounted on 5cms castors with curtain hook curtain springs
- Pre-treated and Epoxy Powder Coated
- Available in Stainless Steel also
- Curtain colour changeable as per requirement.





**CONSOPHARMAPLUS** is an emerging manufacturer-exporter based in New Delhi, capital city of India. The company is founded by highly experienced professionals from Medical Devices and Healthcare industry.

The company has a vision to provide high quality Orthopedic solutions and other healthcare products at a very competitive cost so that the each and every class of society gets quality treatment and care. We ensure our customers and users of our products for the best value of their money.

Our main products include Orthopedic Implants, Orthopedic Instruments, Orthopedic Sets, Surgical Equipments and Rehabilitation Aids etc. The company has production facilities in Gujarat, India which is well known for its industrial expertise. We have completely automated and state-of-art manufacturing unit for orthopedic implants and instruments which assures quality control as per international standards.

Apart from Orthopedic Solutions, we offer our customers a complete solution for their Medical, Surgical and Healthcare needs. These product line include Hospital Furniture, Medical Disposables, Laboratory Supplies, Autoclaves and sterilizers, Suction Units, other medical devices and healthcare equipments.

We welcome dealers, doctors, healthcare professionals to experience all time best services and co-operation. Come, let us work together for serving the community called "Humanity" – a step towards making our world more healthy and beautiful !



CONSOPHARMA PLUS

RZ-C 50, Nihal Vihar, Nangloi, New Delhi - 110041 India.

Call : +91.935050 4074, 987300 4074

Email: [consopharmaplus@gmail.com](mailto:consopharmaplus@gmail.com)